



Barb-to-Female Adapter

Features:



- Standard threads and curvature, no gap, no wobble and good sealing with the corresponding threads
- Imported raw material, non-toxic, tasteless, heat and corrosion resistant.
- PP material is widely used in biomedical, environmental protection and food industries; PTFE material is also suitable for chemical industry and other fields

★ Application Note



Unit: mm

※Can be used with RUNZE flangeless fitting, flanged fitting and barb-to-male adapter to allow for any connection between rigid tube and rigid tube, rigid tube and hose, hose and hose of the same or different tube diameters

Model No.	Tubing ID	Thread	Dimension (mm)			Material
Elbow Barb-to-Female Adapter						
NL05-32-28P	2.4~3.2 mm	1/4-28UNF	A=Φ2.50	B=Φ3.50	C=Φ4.50	PP natural
NL05-32-28E	2.4~3.2 mm	1/4-28UNF	A=Φ2.50	B=Φ3.50	C=Φ4.50	PTFE natural
Straight Barb-to-Female Adapter						
NI06-16-28P	1.0~1.6mm	1/4-28UNF	A=Φ0.99	B=Φ1.60	C=Φ2.90	PP natural
NI06-24-28P	1.6~2.4mm	1/4-28UNF	A=Φ1.59	B=Φ2.40	C=Φ3.80	PP natural
NI06-32-28P	2.4~3.2mm	1/4-28UNF	A=Φ1.97	B=Φ3.50	C=Φ4.50	PP natural
NI06-16-28E	1.0~1.6mm	1/4-28UNF	A=Φ0.99	B=Φ1.60	C=Φ2.90	PTFE natural
NI06-24-28E	1.6~2.4mm	1/4-28UNF	A=Φ1.59	B=Φ2.40	C=Φ3.80	PTFE natural
NI06-32-28E	2.4~3.2mm	1/4-28UNF	A=Φ1.97	B=Φ3.50	C=Φ4.50	PTFE natural